



# Simple answers to complex questions

## Philips Respironics NM3 accessories guide


With a complete selection of adapters and sensors, the Respironics NM3 allows you to monitor adult, paediatric, and neonatal patients. Philips accessories support flexible options for both intubated and non-intubated patients in a variety of settings.

### Key advantages

- In all intubated patients, the Respironics NM3 monitors physiologic gas exchange, deadspace, and alveolar tidal volume with the use of combined flow CO<sub>2</sub> sensors
- In intubated adult patients (VT 200-3000 ml), the Respironics NM3 measures cardiac output
- In non-intubated patients, the Respironics NM3 monitors ETCO<sub>2</sub> and SpO<sub>2</sub> with both mainstream and sidestream capabilities
- VentAssist offers real-time guidance and support for clinical decisions in adult patients

**PHILIPS**  
sense and simplicity

# Capnography accessories for NM3 intubated applications



Capnostat 5  
mainstream  
CO<sub>2</sub> sensor

Single-patient  
combination sensors  
(Neonatal)

Single-patient  
airway adapters  
(Neonatal)

Cardiac  
output sensor  
(Standard)

## Capnostat 5 mainstream for long-term use

Description	Part Number
Capnostat 5 mainstream CO <sub>2</sub> sensor	1015928
<b>Single-patient combination CO<sub>2</sub> flow sensors</b>	
Adult	9767
Paediatric	976600
Neonatal	976500
<b>Single-patient airway adapters</b>	
Adult/paediatric	6063
Neonatal	6312
<b>Cardiac output sensor</b>	
Small (200-500 ml)	895000
Standard (400-1000 ml)	8951
Large (750-1500 ml)	895200



© 2009 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

Printed in Australia  
4522 962 56101 \* OCT 2009

# Teleflex

Proudly Distributed by **Teleflex Medical**

**Teleflex Medical Australia Pty Ltd**  
Level 4 Building B · 201 Coward Street  
Mascot 2020 NSW · Australia  
Ph: 1300 360 226 · [austCS@teleflex.com](mailto:austCS@teleflex.com)  
[www.teleflexmedical.com.au](http://www.teleflexmedical.com.au)

**Teleflex Medical New Zealand**  
12 Victoria Street  
Lower Hutt 50103 · New Zealand  
Ph: 0800 601 100 · [nzsupport@teleflex.com](mailto:nzsupport@teleflex.com)  
[www.teleflexmedical.co.nz](http://www.teleflexmedical.co.nz)